# Flowtron

## DVT40

ARJOHUNTLEIGH  
GETINGE GROUP

**EN** To be used only under the direction of a physician - Non-Sterile -  
Not made with natural rubber latex - For Single Patient Use Only

**Caution:** Federal law restricts this device to sale by or on the order of a licensed practitioner.

**ES** *For uso sólo con Flottron® DVT-prevention pumps manufactured by ArjoHuntleigh. Must not be used with Flottron® Hydroven 3 o Flottron® Hydroven 12 Intermittent Pneumatic Compression (IPC) pumps.*

**DE** *More comprehensive information on the garment can be found in the relevant Flowtron DVT pump Instructions For Use document.*

### Instructions for Use

1. Plug the pump into a suitable electrical outlet. Do not turn the pump on at this time.
2. Remove the garments from the sealed bag. The garments may be used on either leg. Unfold the garment and position the inflatable bladder along the leg.
3. Snugly wrap the garment around the patient's leg and secure the fastener tabs. Repeat for the other leg.
4. Attach the garment to the pump tubing set ensuring a 'click' is heard from each snap-lock connector.
5. If your Flowtron DVT pump has an adjustable pressure regulator, turn the regulator dial to the recommended pressure unless otherwise directed by the physician. For further details, refer to the relevant Flowtron DVT pump Instructions For Use document.
6. Turn the pump on. The green power indicator lights should illuminate.

**Indications**  
Intermittent Pneumatic Compression (IPC) is indicated for use for the prevention of deep vein thrombosis (DVT) due to the presence of risk factors for thrombus formation during orthopaedic, trauma, urologic and general surgery. It is also suitable for other patient groups including neurology, critical care, general medical and obstetrics.

### Recommendations

General recommendations:

- Garments should be removed regularly to inspect the skin
- Where appropriate, patients should be instructed in the proper use of the system, the purpose of therapy and any problems that should be reported to the nursing staff

**DE** *® and ™ are trademarks belonging to the ArjoHuntleigh group of companies  
As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.*

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© og ™ er varemærker tilhørende ArjoHuntleigh-koncernen